


METHOD FEEDER LEAD SYSTEM


This unique lead/feeder system can be used for making 'The Method' with the Wolf Advance Ball Maker. The plastic feeder is in two halves, one half is filled with lead/steel discs and felt discs (used to soak up attractor liquid) and is clipped together with the swivel clipped in to the top position. The feeder which has small holes for leaking out the liquid attractor slowly can be injected with the liquid before screwing on to the Ball Maker. The Ball maker is then inserted in to the prepared method mix and is compressed to make the aerodynamic egg shaped method balls wrapped around the feeder lead. The Method/Feeder is then attached to the lead clip and because of the aerodynamic shape and the weight of approx. 50/60 grams it can be cast long distances.



Attach swivel to peg on Feeder.




Insert Ball Maker into prepared Method Mix.



Discs and foam inserts provided. Spare leads and foam discs can be bought separately.



Compress Ball Maker to make the Egg-shaped Method Ball wrapped around the Method Feeder.



Depending what weight required different combinations of lead and foam inserts can be used.



Clip two halves carefully together.



Unscrew and attach to safety lead clip. Embed hook in to Method Ball to avoid tangles.



Liquid Attractor can be squeezed or injected into feeder.




Liquid attractor leaking.



Liquid Attractor will leak out slowly as the method mix breaks down.



Screw Method Feeder on to Ball Maker.



RESULT!